



Audiometric Data Management

MINIMUM SYSTEM REQUIREMENTS

Minimum PC requirements

- CPU: 1.6 GHz
- RAM: 1 GB
- Hard Disk: 20 MB
- USB Ports: 2
- Minimum Screen Resolution: 1024 x 768

Supported Operating Systems Editions

- Windows® 7 Professional
- Windows® 8 Pro
- Windows® 10 Pro

Operating Systems Languages

- English
- Simplified Chinese
- Russian
- Italian
- Japanese
- Portuguese
- French
- German
- Spanish
- Korean

Other Required Applications (included)

- .NET 4.5
- Adobe Acrobat reader software
- GSI Instrument drivers

Optional Applications

- OtoAccess™
- Noah® 4